

Climate chambers / HUMIDITY CHAMBERS HCP

Options	Order No.	50	105	150	240
Voltage 115 V, 50/60 Hz	X2			●	
Battery-buffered ControlCOCKPIT: uninterrupted supply for the entire display unit (ControlCOCKPIT) and therefore complete documentation of all parameters even when there is a power failure	C2			●	
Entry port with 23 mm clear diameter, can be closed, standard position	left centre/top			●	
	right centre/top			●	
4 - 20 mA current loop interface	Temperature controller, actual value (0 to +100 °C = 4 - 20 mA)			●	
	Humidity controller, actual value (0 to 100 % rh = 4 - 20 mA)			●	
Door with lock and key (safety lock)	B6			●	
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)	H6			●	
MobileALERT, notification by SMS in case of any error or alarm of the device (requires option H6)	C3			●	
MobileALERT for 2 alarm notifications; notification by SMS. temperature and humidity alarm (requires option H6)	C4			●	

Accessories	50	105	150	240
Stainless steel grid, electropolished	E20164	E20165		E43118
Perforated stainless steel shelf	E35160	E37418		E35158
Central water supply with filter cartridges for connection to the domestic water supply. Product information on demand		ZWVR6		
Central water supply without filter cartridges for connection to the domestic water supply (only for demineralised water with a conductivity of 5 to 10 µS/cm and a pH value between 5 and 7). Product information on demand		ZWVR7		
Stacking set (consisting of stacking corners, one connecting plate for the rear, two wall brackets) for stacking of two units of same size (requires subframe B33509 or B43598)		-		B48129
Stacking set (consisting of stacking corners, one connecting plate for the rear, two wall brackets) for stacking of two units of same size			B42114	
Stacking set (4 pcs) for stacking of appliances of same size	B29744			
Subframe, on castors (height 120 mm; stainless steel, material 1.4301)				B43598
Subframe (622 mm high) adjustable in height (sizes 150/240: should not be used for 2 stacked units)	B33504	B33505		B33506
Subframe (130 mm high); sizes 150/240: only in combination with the corresponding stacking sets for stacked appliances	B33507	B33508		B33509
External measuring instrument with additional measuring head for temperature and humidity measurement. Product information on demand		B04714		

Accessories Software and Connectivity	50	105	150	240
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Respective IQ/OQ documents available in German and English language (without surcharge)			FDAQ1	
Integration of additional units (up to max. 31 units) into an already existent FDA-software licence			FDAQ2	
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number			B33170	
Ethernet connection cable 5 m for computer interface			E06189	
USB-Ethernet adapter			E06192	