

Medical devices / INCUBATORS Im / models INm, INmplus, IFm, IFmplus

Options	Order No.	30	55	75	110	160	260	450	750
Voltage 115 V, 50/60 Hz	X2					●			
Extended overtemperature protection by additionally integrated Pt100 sensor for independent temperature monitoring for models with SingleDISPLAY	A6					●			
Temperature restriction; Temperatures: +50, +55, +60, +65, +70 or +75 °C (Please indicate upon ordering)	A8					●			
Accessories	30	55	75	110	160	260	450	750	
Stainless steel grid, electropolished	E28884		E20164		E20165	E28891		E20182	
Reinforced stainless steel grid, electropolished, max. loading 60 kg. Please consider max. loading of chamber					E29767	E29766			
Perforated stainless steel shelf	B29727		B03916		B00325	B29725		B00328	
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution)	E02070		E02072		E02073	E29726		E02075	
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution)	B04356		B04358		B04359	B29722		B04362	
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), straight, for exhaust air ducting (if necessary for connection by hose)					B29718				
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), angled, for exhaust air ducting (if necessary for connection by hose)					B29719				
Wall bracket for wall mounting	B29755	B29756	B29757	B29758	B29759				
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29728	B29730	B29732	B29734	B29736	B29738	B29740	B29742	
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29729	B29731	B29733	B29735	B29737	B29739	B29741	B29743	
Stacking set (4 pcs) for stacking of appliances of same size			B29744						
Set of height adjustable feet (4 pcs)			B29768						
Subframe, adjustable in height (size 30 to 75: height 600 mm, from size 110: height 500 mm)	B29745		B29747		B29749	B29751	B29753		
Subframe, on castors (size 30 to 75: height 660 mm, size 110 to 160: height 560 mm)	B29746		B29748		B29750				
Subframe, adjustable in height, height 130 mm	B33657		B33659		B33661	B33664			
Accessories Software and Connectivity	30	55	75	110	160	260	450	750	
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit (only for units with TwinDISPLAY). Respective IQ/OQ documents available in German and English language (without surcharge)					FDAQ1				
Integration of one additional unit (up to max. 31 units) into an already existent FDA-software licence (only for units with TwinDISPLAY)					FDAQ2				
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number (only for units with TwinDISPLAY)					B33170				
USB stick with documentation software AtmoCONTROL and operation manual for products with SingleDISPLAY (the standard equipment of appliances with TwinDISPLAY includes one USB stick with AtmoCONTROL). When reordering please specify serial number					B33172				
Ethernet connection cable 5 m for computer interface					E06189				
USB-Ethernet adapter					E06192				

Standard equipment included in the standard scope of delivery. Prices for options are only valid when ordering new appliances.
Not all options/accessories are combinable with each other. Please contact us for individual combination requests.