

Medical devices / UNIVERSAL OVENS Um / models UNm, UNmplus, UFm, UFmplus

Options	Order No.	30	55	75	110	160	260	450	750	1060
Voltage 115 V, 50/60 Hz	X2				●					
Extended overtemperature protection by additionally integrated Pt100 sensor for independent temperature monitoring for models with SingleDISPLAY	A6					●				
Temperature restriction; Temperatures: +60, +70, +80, +95, +100, +120, +160, +180, +200 or +250 °C (Please indicate upon ordering)	A8					●				

Accessories	30	55	75	110	160	260	450	750	1060
Stainless steel grid, electropolished	E28884	E20164		E20165		E28891	E20182		B41251
Reinforced stainless steel grid, electropolished, max. loading 60 kg. Please consider max. loading of chamber				E29767		E29766			B32550
Perforated stainless steel shelf	B29727	B03916		B00325		B29725	B00328		
Reinforced stainless steel shelf, max. loading 60 kg. Please consider max. loading of chamber									B32549
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution)	E02070	E02072		E02073		E29726	E02075		B32599
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution)	B04356	B04358		B04359		B29722	B04362		B29769
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), straight, for exhaust air ducting (if necessary for connection by hose)						B29718			
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), angled, for exhaust air ducting (if necessary for connection by hose)						B29719			
Wall bracket for wall mounting	B29755	B29756	B29757	B29758	B29759				
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29728	B29730	B29732	B29734	B29736	B29738	B29740		B29742
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29729	B29731	B29733	B29735	B29737	B29739	B29741		B29743
Stacking set (4 pcs) for stacking of appliances of same size		B29744							
Set of height adjustable feet (4 pcs)				B29768					
Subframe, adjustable in height (size 30 to 75: height 600 mm, from size 110: height 500 mm)	B29745	B29747		B29749		B29751	B29753		
Subframe, on castors (size 30 to 75: height 660 mm, size 110 to 160: height 560 mm)	B29746	B29748		B29750					
Subframe, adjustable in height, height 130 mm	B33657	B33659		B33661		B33664			

Accessories   Software and Connectivity	30	55	75	110	160	260	450	750	1060
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit (only for units with TwinDISPLAY). Respective IQ/OQ documents available in German and English language (without surcharge)						FDAQ1			
Integration of one additional unit (up to max. 31 units) into an already existent FDA-software licence (only for units with TwinDISPLAY)						FDAQ2			
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number (only for units with TwinDISPLAY)						B33170			
USB stick with documentation software AtmoCONTROL and operation manual for products with SingleDISPLAY (the standard equipment of appliances with TwinDISPLAY includes one USB stick with AtmoCONTROL). When reordering please specify serial number						B33172			
Ethernet connection cable 5 m for computer interface						E06189			
USB-Ethernet adapter						E06192			

Standard equipment included in the standard scope of delivery. Prices for options are only valid when ordering new appliances. Not all options/accessories are combinable with each other. Please contact us for individual combination requests.