



COLONY COUNTER

Code BRJ003, BRJ004

The Nahita-Blue digital colony counter has been developed to ensure consistent results and reliable operation in routine microbiological counting tasks, facilitating proper management and recording of the data obtained.

■ **Pressure sensor design**

Allows counting using any marker pen, preventing deterioration or contamination of the sample.

■ **Quadruple display panel**

Displays up to four counts simultaneously: the current count and the three previous counts.

■ **Advanced counting functions**

1. Automatic averaging: calculates the average of multiple counts with a single press.
2. Countdown function: allows correction of excessive counts made by mistake.
3. Quick reset function.

■ **Compatibility**

Suitable for Petri dishes of different diameters: 60 mm (with adapter ring) and 90–150 mm.

■ **Optional software**

Facilitates data traceability with functions such as delete, search, print, and save data in Excel (.csv) format.



Technical specifications

Code	BRJ003
Display	LED
Magnification	2x
Operating functions	Average, zero reset, count back
Illumination	LED (adjustable)
Pressure sensitivity	Adjustable
Acoustic signal	Adjustable
Output	USB
Weight	5.6 kg
Dimensions	35x30x37 cm

Optional software

Code	BRJ004
Description	Software for digital colony counter