



► PARAFFIN DISPENSER - ZFP013

- 1 | Easy operation through a 6-button control panel.
- 2 | Power LED indicator.
- 3 | Backlight LCD screen.
- 4 | Heating tank with lid and a capacity of 5 L of paraffin.
- 5 | Maximum heating temperature: 99° C.
- 6 | Dispenser tap operated by a button.
- 7 | Possibility of programming heating start and end time.



| | |
|-------------------------|-----------------------------------|
| Reference | ZFP013 |
| Deposit capacity | 5 L |
| Programmable parameters | Start time and end of temperature |
| Dimensions | 410x440x360 mm |
| Electrical power supply | 220 V, 50 Hz |