

# Universal suction kit

## BioDolphin with different adaptors

### Ideal for:

- cell culture aspiration
- laboratory waste suction

### Features:


- user-friendly
- continuous suction without pressing the button
- quick connection between fitting and handle
- storage rack to keep order
- built-in liquid collecting tray
- autoclavable (except storage rack)



BioDolphin


Model	BioDolphin complete kit
Operation	Press to suck
Port thread	3/16 inch (5 mm)
Dimension (LxWxH)	15 x 15.2 x 14.5
Order number	00829 199100-00

Different waste bottles can be found on page 193  
Vacuum pumps can be found from page 190 to 196



# Universal suction kit

Complete kit contains following material:



#	Model	Material	separate order number
1	Handle grip	POM	00829 199100-30
2	Rack		00829 199100-61
3	1-Channel adaptor (80 mm)	Stainless steel, POM	00829 199100-48
4	1-Channel adaptor (40 mm)	Stainless steel, POM	00829 199100-47
5	1-Channel tip adaptor	POM	00829 199100-46
6	8-Channel adaptor (40 mm)	Stainless steel, POM	00829 199100-49
7	8-Channel tip adaptor with ejector	POM	00829 199100-50
8	1-Channel tip adaptor with ejector	POM	00829 199100-75

