

## DRYSC43 Laboratory Drying Oven Economic 43L

### Product Details

Manufacturer: Colo Lab Experts

**COLO CODE: DRYSC43**



### Technical Specification

- Stainless-steel chamber
- PID controller
- Uniform distribution of air temperature
- Forced-air convection
- Double Layer glass door, larger viewing window
- Temperature range RT-300°C
- Temperature control: 1°C
- Temperature Stability:  $\pm 1^\circ\text{C}$
- Timer: 1-9999min
- Shelves: 2pcs

Model	DRYSC43C	DRYSC71C	DRYSC136C	DRYSC225C
Power Supply		220V/50Hz		
Temperature range		RT+5 - 300°C		
Chamber size mm	350*350*350	450*450*350	550*550*450	760*600*500
Power Kw	1,2	1,6	2	5
Product size	660*530*580	760*530*580	860*530*780	870*670*1030
Packing size	760*620*660	820*620*760	960*720*860	1030*720*1120
Net weight	44	54	58	110
Gross weight	56	66	79	140

### Optional Accessories

#### Optional accessories:

Shelves, RS485 interface, printer, recorder, external communication, remote control, program temperature control device, wireless SMS alarm, USB data storage.