

DRYSD30 Laboratory Drying Oven Standard 30L

Product Details

Manufacturer: Colo Lab Experts

COLO CODE: DRYSD30



Technical Specification

Well coordinated bottom air fan (fan temperature less than 50°C) and heater and vertical hot air circulation ensure the temperature evenness. Stainless steel inner chamber & removable shelves. PID intelligent double digital temperature controller with functions of timing and over-temperature protection. Door with wide double-glass observation window, adjustable airtight buckle lock which efficiently prevent heat loss. Low noise fan with the air channel keeps good temperature uniformity.

Volume: 30 L
 Inner Chamber Size: 310 x 310 x 310 mm
 Voltage: 220V / 110V
 Power: 0.8 KW
 Temperature Motion: +1 C
 Temperature Range: RT+10 C to 300 C
 Overall Size: 550 x 575 x 780 mm
 Timer: 0 - 9999 min
 Over heat protection

Model	DRYSD30	DRYSD45	DRYSD65	DRYSD125
Power Supply	220V/50Hz			
Temperature range	RT+5 - 300°C			
Chamber size mm	310*310*310	350*350*350	400*360*450	500*450*550
Power Kw	0,8	1,2	1,6	2,3
Product size	450*490*685	490*530*730	540*540*830	650*630*930
Packing size	550*570*580	630*600*820	650*640*920	750*730*1020
Net weight	33	37	44	60
Gross weight	40	45	54	72

Optional Accessories