


MAIN FEATURES





- > Adjusted for use in cytology, urology, microbiology, virology, oncology
- > Wide range of consumables for different applications
- > Centrifugation of up to 12 samples
- > Robust design
- > Microprocessor controlled functions of cover, time and speed
- > Adjustable acceleration and braking system
- > Special START - STOP features
- > FAST SPIN function
- > Possibility of programming



CYTOLOGY SOLUTIONS	CENTRIC 250	CENTRIC 350
No. of places in rotor	8	12
Max. speed	2.500 rpm	2.500 rpm
Max. centrifugal force	630 x g	860 x g
Brushless motor	YES	YES
Noise level (dBA)	48	49



ACCESSORIES FOR CytoCENTRIC CENTRIFUGES

ROTOCLIPS STAINLESS STEEL JC 390	ROTOFUNELL Reusable JC 406	FILTER CARDS JC 407	SLIDES JC 311	ECOfunnel® Reusable JC 313	ECOseal JC 315	ECOslides JC 314	Filter cards 345-100
Order ref. No.: 312491	Order ref. No.: 496090	Order ref. No.: 484130	Order ref. No.: 484136	Order ref. No.: 334724	Order ref. No.: 334729	Order ref. No.: 334733	Order ref. No.:
4 pcs in a pack	12 pcs in a pack	200 pcs in a pack	100 pcs in a pack	4 pcs in a pack	100 pcs in a pack	100 pcs in pack	100 pcs in pack
	 For sample volumes of up to 0.5ml Sedimentation field with a diameter of 6mm		 For preparing large-volume samples of up to 6ml Sedimentation field of 22x15mm for a higher cell yield Only for rotor with 12 places				

We reserve the right to alter specification details etc. without prior notice or liability!

V:05-23