

COOL INCUBATOR 56 L

INCUBATORS

A digital benchtop incubator with both cooling and heating capabilities. Experience outstanding temperature stability, low energy consumption, and a compact design. Perfect for microbiology and hematology applications, this incubator ensures your experiments are in the ideal environment.

TECHNICAL SPECIFICATIONS

Model	Cool Incubator 56 L
Temperature set range	ambient -10°C to 50°C
Controls	digital
Stability	±0,2°C at 37°C
Accuracy	±0,5°C
Capacity	56 liters
Shelves	5 positions, 2 shelves included
Display	LED
Max. power consumption	160 W
Chamber Dimension (H x W x D)	400 x 394 x 354 mm
Dimension (H x W x D)	580 x 480 x 525 mm
Crating	600 x 500 x 650 mm
Weight (Net)	35 kg
Permitted working ambient	temperature 2°C to 35°C, rh 75% non-condensing

ORDER INFORMATION

230 V, 50/60 Hz	Cool Incubator 56
Code	225.000
Order ref.	528457
120 V, 50/60 Hz	Cool Incubator 56 120 V
Code	1.225.000
Order ref.	534501

Additional shelf (for 56 liter incubator)

Code	781.1.560
Order ref.	377193
Shelf autoclavable	yes

