



BASIC line

Article Number: 23-201267-030

Product Range Safety storage cabinet
Product Category Type 30 - flammable liquids

General

Line	BASIC
Version	standard
Size	XL
Type	30
Storage type	passive storage
Functionality	Storage
Earthing connection	Yes
Connection socket	Yes

Test criteria and documentation

DIN EN 14470-1	Yes
DIN EN 16121/16122	Yes
CE conformity	Yes

Technical specifications

Outer width	1194 mm
Outer depth	612 mm
Outer height	2045 mm
Max. load bearing capacity per storage area	75 kg
Max. load	360 kg
Max. total weight	810 kg
Max. volume of single container	30 L
Collection volume	33 L
Weight (product)	450 kg

Cabinet

Material of cabinet cladding	Sheet steel
Cabinet colour	RAL 7035
Material of inner lining	Decor panels
Colour of inner lining	RAL 7035
Base	Classical
Base height	89 mm
Base colour	RAL 7015

Doors

No. of doors	2
Kind of door	Wing door
Material of door cover	Sheet steel
Door colour	RAL 7035
Door operation	Classical
Door technology	Classical
Door handle	Bow handle matt nickel
Lock	Cylinder lock
Locking system	Open
Locking points	1-point

Interior fittings

Material of interior fittings	Sheet steel
Colour of interior fittings	RAL 7035
Bottom tray	Fixed
Storage shelf	Yes
Storage levels	3
Storage areas	3
setting height	526 mm