

## POCKET'S

EB02.036

### Molle side pockets

This pair of auxiliary pockets can be attached to any item with a molle system, such as the paramedical rescue kit PARAMED'S, giving extra capacity for auxiliary material.



### Description

- Made of 1000D washable polyamide, they have a bellows-type rack that can be opened in case you need to expand their capacity.
- Its front has molle system which allows them to anchor to them other compartments.
- Wide handles with reflective night-day band for easy viewing in an emergency situation.

### Additional information

Weight	0.18 Kg
Dimensions	6 x 16 x 28 cm
Material	1000D POLYESTER
Color	RED
Capacidad (Litros)	2.39
Contenido máximo recomendado (Kg)	4.00