

GEL

EB09.007

Reusable cold gel

Perfect size gel to include in our ampoule holders. Thought to be able to keep the medication cool at scheduled departures or emergencies. Indications: Haematomas, sprains, injuries, high temperature and whenever cold application is required. Put the gelatine bag in the freezer for 6/12h. Take it out, wrap it up in a cloth and apply where it requires.



Features

Applications: Bruising, sprains, hyperthermia and in all cases where cryotherapy is necessary.

Technical data

Weight	0.12 Kg
Dimensions	1 x 17 x 12 cm
Material	GEL
Color	BLUE
Capacity (Liters)	0.19