

EXTREME'S

EB02.008

Emergency bag. Basic life support

It is a basic life support (BLS) bag. Its honeycomb internal structure allows the user to easily organize all the contents. It has several detachable compartments in different colors.



© Marcas, Diseños y presentación comercial protegidos ELITE BAGS S.L., 2018
 © ALL TRADE MARKS, DESIGNS & TRADE DRES S. ELITE BAGS S.L., 2018

Description

- Padded handle
- Reflective night/day high visibility bands
- Reflective parachutist rope handles that glow in the dark
- Ring for descents in air transport
- Oxygen tube extraction washer (for disposable oxygen bottles)
- Cardholder
- Removable backpack straps
- Two large capacity external pockets with elastic bands and pockets with transparent window
- Fully adaptable interior thanks to its removable velcro padded dividers
- Reinforced base
- Four removable compartments with transparent front with the international color code
- Isothermal ampoule with capacity for 27 ampoules
- Interior adaptable thanks to its separators with velcro system
- Holding straps and space for the secretion aspirator
- Velcro identifiers in the side pockets
- Velcro identifier with the Cross of Life

[embed]https://www.youtube.com/watch?v=UoXYzzPHtTM[/embed]

Additional information

Weight	2.86 Kg
Dimensions	28 x 47 x 22 cm
Material	1000D POLYESTER
Color	RED
Contenido máximo recomendado (Kg)	12.00
Capacidad (Litros)	28.95L