

MEDIC'S EB06.011

Sports medical bag



It is a sports medical bag with a practical and comfortable design. Its honeycomb internal lay out allows to quickly find the required material at any sports or medical emergency situation.

Features

1. Padded handle with velcro fastener.
2. Cardholder.
3. Reinforced structure thanks to a webbing strap.

Technical data

- **Measures:** 50 x 25 x 31 cm
- **Capacity:** 38,75 L
- **Weight:** 2.23 kg
- **Material:** 1200D Polyester. Washable
- **Colour:** Black
- **Maximum recommended contents weight:** 12 kg