

# ELECTROMAGNETIC STIRRERS

## Stirrers without heating



**F20**  
SMALL VOLUMES



Useful in the microbiological and biochemical fields



The stainless steel plate facilitates cleaning and sterilization



Electronic speed control DC motor

### Technical data

### F20

Stainless steel working surface LxD	mm	110x110
Max stirring capacity with H <sub>2</sub> O	lt	2
Variable speed	Rpm	100-1100
Hole for rod	M 8	No
Motor power	W	10
Dimensions LxWxH	mm	110x115x60
Weight	kg	0,7
Stirring bar suggested Ø	mm	6x20

### Code

**601.0101.20**

V/HZ 220/240-50/60

SAFETY CLASS 0

PROTECTION CLASS IP 41