



Mini-incubators of 28 lt – ICT 28

Summary

28 lt mini-incubators Falc are small instruments for microbiologic applications. They are used to monitor the temperature stability of samples at 37°C using natural ventilation not to compromise them.

Differently from big volume incubators, mini-incubators have plexiglass door to check samples without opening it and their small sizes make them be easy-friendly products, spread in all labs.

Working temperature is below 80°C.

Falc incubators create the correct microclimate suitable for microorganisms' cultures without changing surrounding conditions.

Description

*Natural ventilation

*Plexiglas door

*Digital temperature control

* 10 mm hole placed on the upper part for insertion of probe or thermometer



Features

Item	720.2070.14
Temperature	80 °C
Precision	± 0.5
Risolution	0.1
Heating power	W 420
Internal dimension (LxWxH)	mm 400x250x300



Number racks included	n. 2
Max racks number	n. 5
External dimension (LxWxH)	640x360x440 mm
Power	V/Hz 230-50/60
Weight	Kg 26.5
Security class	0
Protection class	IP 42
Capacity	lt 28.0