



Basic Line- Natural ventilation oven 120 lt - STE-N 120

Summary

Natural ventilation basic oven of 120 lt

Basic ovens are reliable and robust as standard Falc models (STF-N/STF-F) but they are made with simple structural solutions to get them cheaper.

Falc Basic Ovens are perfect instruments for generic thermostatic purposes such as glassware drying and sterilizing that don't need specific temperature. To do this, they are used with or without ventilation. Unlike incubators, ovens reach higher temperatures, up to 300°C.

Description

- * Natural ventilation

- * External steel structure varnished with epoxy paint
- * Internal chamber in 18/8 stainless steel

- * Security class 0
- * 45 mm hole placed on the posterior part for fume discharge, with regulation device for opening/closing passage
- * 10 mm hole placed on the upper part for insertion of probe or thermometer

- * Manual chimney opening



Features

Item	720.2079.05
Min temperature	Ambiente +5 °C
Max temperature	+300 °C
Precision	150°C ±3,0 / 300°C ±3,5
Risolution	0,1°C
Heating power	W 1600
Internal dimension (LxWxH)	390x350x390 mm
Number racks included	n. 2
External dimension (LxWxH)	820x710x740 mm
Power	V/Hz 220/240 - 50/60
Power	V/HZ 220/240-50/60
Weight	Kg 68



Maximum number racks	n° 6
Security class	0
Protection class	IP 42
Capacity	lt 120