



---

## Mini analog oven 28 lt - STZ 28

### Summary

Mini analog oven of 28 lt

Natural convection drying ovens with notable thermal stability.

Specially indicated for laboratories and pharmacies that need to sterilize small equipment.

### Description

\*Natural convection drying ovens with notable thermal stability. Specially indicated for laboratories and pharmacies that need to sterilize small equipment.

\*Internal chamber in aluminum for a better heat diffusion.

\*Frontal panel with bimetallic fluid expansion probe thermoregulator

\*Analogue timer (mechanical) from 0-120 min. or infinite

\*Natural ventilation

\* 10 mm hole placed on the upper part for insertion of probe or thermometer



## Features

Item	720.2070.08
Temperature	150 °C
Precision	±5
Heating power	W 600
Internal dimension (LxWxH)	mm 400x250x300
Number racks included	n. 2
Max racks number	n. 5
External dimension (LxWxH)	640x360x440 mm
Power	V/Hz 230-50/60
Weight	Kg 26.5
Timer	0-120 or (∞)
Protection class	IP 40
Capacity	lt 28.0