



Mini analog oven 5.4 lt - STZ 5.4

Summary

Mini analog oven of 5.4 lt

Natural convection drying ovens with notable thermal stability.
Specially indicated for laboratories and pharmacies that need to sterilize small equipment.

Description

*Natural convection drying ovens with notable thermal stability. Specially indicated for laboratories and pharmacies that need to sterilize small equipment.

*Internal chamber in aluminum for a better heat diffusion.

*Frontal panel with bimetallic fluid expansion probe thermoregulator

*Analogue timer (mechanical) from 0-120 min. or infinite

*Natural ventilation

* 10 mm hole placed on the upper part for insertion of probe or thermometer



Features

Item	720.2070.04
Temperature	150 °C
Precision	± 5
Heating power	W 360
Internal dimension (LxWxH)	mm 190x150x190
Number racks included	n. 2
Max racks number	n. 4
External dimension (LxWxH)	400x260x300 mm
Power	V/Hz 230-50/60
Weight	Kg 10
Timer	0-120 or (∞)
Protection class	IP 40
Capacity	lt 5.4