



Mini digital oven 28 lt - STZ-D 28

Summary

Mini Digital Oven of 28lt.

Natural convection drying ovens with notable thermal stability.

Specially indicated for laboratories and pharmacies that need to sterilize small equipment.

Description

*Natural convection drying ovens with notable thermal stability. Specially indicated for laboratories and pharmacies that need to sterilize small equipment.

*Internal chamber in aluminum for a better heat diffusion.

*Frontal panel with bimetallic fluid expansion probe thermoregulator

*Digital display

*Natural ventilation

* 10 mm hole placed on the upper part for insertion of probe or thermometer



Features

Temperature	+5-amb - +200°C
Precision	+/-2 °C
Risolution	0.1
Heating power	600 W
Internal dimension (LxWxH)	400x250x300mm
Number racks included	2
Max racks number	5
External dimension (LxWxH)	640x360x440 mm
Weight	26.5 kg
Capacity	28 Lt