



Mini digital oven 5.4 lt - STZ-D 5.4

Summary

Mini digital oven of 5.4 lt

Natural convection drying ovens with notable thermal stability.

Specially indicated for laboratories and pharmacies that need to sterilize small equipment.

Description

*Natural convection drying ovens with notable thermal stability. Specially indicated for laboratories and pharmacies that need to sterilize small equipment.

*Internal chamber in aluminum for a better heat diffusion.

*Frontal panel with bimetallic fluid expansion probe thermoregulator

*Digital display

*Natural ventilation

* 10 mm hole placed on the upper part for insertion of probe or thermometer



Features

Temperature	+5 amb - 200°C
Precision	+/-2°
Risolution	0.1
Heating power	360
Internal dimension (LxWxH)	190x150x190mm
Number racks included	2
Max racks number	4
External dimension (LxWxH)	400x260x300 mm
Weight	10 kg
Capacity	5.4 lt