










LAUDA Varioshake accessories

For every purpose

	Designation	VS 8 OE / BE / O / B
Shaking trays, supports, rack attachments (suitable for/max. number)		
	Adhesive mats, black, 200 x 200 mm, can be cut to size	2
	Non-slip support	1
	Non-slip support	-
	Shaking tray	1
	Shaking tray	-
	Shaking tray	-
	Shaking tray	-
	Universal attachment	1
	Universal attachment	-
	Universal attachment	-
	Platform frame	-
Clamp for screwing onto trays (suitable for/max. number)		
A000044		
	Clamp for Erlenmeyer flask 25 ml	45
	Clamp for Erlenmeyer flask 50 ml	25
	Clamp for Erlenmeyer flask 100 ml	16
	Clamp for Erlenmeyer flask 200 ml	12
	Clamp for Erlenmeyer flask 250 - 300 ml	9
	Clamp for Erlenmeyer flask 500 ml	9
	Clamp for Erlenmeyer flask 1000 ml	4
	Clamp for Erlenmeyer flask 2000 ml	2
	Clamp for 50 ml separating funnel (Squibb)	6
	Clamp for 100 ml separating funnel (Squibb)	6
	Clamp for 250 ml separating funnel (Squibb)	4
	Clamp for 250 ml separating funnel (conical)	4
	Clamp for 500 ml separating funnel (Squibb)	3
Clamp/ mount for screwing onto trays (suitable for/max. number)		
A000044		
	Test tube rack for max. 24 tubes of 12 - 17 mm Ø, length 75 - 160 mm	3
	Test tube rack for max. 16 tubes of 25 - 29 mm Ø, length 75 - 160 mm	2
	Mount for test plates	4

	VS 15 T/ R	VS 15 O/ B	VS 30 O	VS 20 OH	VS 60 OI	VS 45 OI	VS 150 OI	Catalog number
	-	4	8	-	4	6	8*	A000041
	-	-	-	-	-	-	-	A000042
	1	1	-	-	1	-	-	A000043
	-	-	-	-	-	-	-	A000044
	1	1	-	-	1	-	-	A000045
	-	-	-	-	-	2	-	A000046
	-	-	1	-	-	-	2	A000047
	-	-	-	-	-	-	-	A000048
	1	1	-	-	1	-	-	A000049
	-	-	1	-	-	-	-	A000050
	1	-	-	-	-	-	-	A000051
	A000045	A000045	A000047	-	A000045	A000046*	A000047*	
	79	79	99	-	79	52	99	A000025
	49	49	99	-	49	33	99	A000026
	36	36	50	-	36	22	50	A000027
	22	22	26	-	22	15	26	A000028
	16	16	26	-	16	13	26	A000029
	12	12	26	-	12	10	26	A000030
	9	9	12	-	9	6	12	A000031
	4	4	9	-	-	3	9	A000053
	-	11	-	-	-	-	-	A000054
	-	11	-	-	-	-	-	A000055
	-	8	-	-	-	-	-	A000056
	-	8	-	-	-	-	-	A000057
	-	6	-	-	-	-	-	A000058
	-	A000045	A000047	-	A000045	A000046*	A000047*	
	-	6	9	-	6	3	9	A000059
	-	-	-	-	-	-	-	A000060
	-	6	-	-	6	6	15	A000061

* Specified maximum quantity per tray
 VS 45 OI (tray A000046): Only 1 tray can be used from a flask size of 300 ml
 VS 150 OI (tray A000047): Only 1 tray can be used from a flask size of 1000 ml