

Heating bath circulators



SB series

- For external and internal application temperature
- Microprocessor pid control
- Powerful circulation pump to keep uniform temperature of the bath
- Bottom of the tank with rounded corner for easy cleaning


Technical data		SB5	SB15	SB24
Tank Capacity	lt	5	12	18,5
Drain Cock		no	no	yes
Digital timer	hh:mm	∞ o 99:59	∞ o 99:59	∞ o 99:59
Adjustable temperature	°C	+5 amb. ÷ 150	+5 amb. ÷ 150	+5 amb. ÷ 150
Precision temperature	°C	±0,5	±0,5	±0,5
Display resolution	°C	0,1	0,1	0,1
Pump range mbar 190	lt/min	9	9	9
Discharge head mt. 1.8	lt/min	7,8	7,8	7,8
Heating power	W	1300	1300	1300
Dimension on the bottom WxDxH	mm	250x100x150	290x255x150	440x560x200
Tank dimensions WxDxH	mm	160x150x150	190x300x150	360x300x150
External dimension WxDxH	mm	330x180x370	360x350x370	530x350x370
Weight	kg	7	8	11
Item		612.1029.05	612.1029.15	612.1029.24

ACCESSORIES



SERPENTINE SET
To work up to -10 ° C
by an external cooling system

Item
612.1029.01



Stainless steel test - Tube rack with 3 shelves

	1 for tank	1 for tank	2 for tank
	12 holes Ø21	48 holes Ø21	48 holes Ø21
Item	212.1250.02	212.1250.05	212.1250.05
	20 holes Ø17	60 holes Ø17	60 holes Ø17
Item	212.1250.03	212.1250.06	212.1250.06

PUMP SET
It allows the external power
supply via a closed circuit

Item
612.1029.02

COMPLETE SET
Serpentine set + pump set

Item
612.1029.03