



ML0423WN



ML0423WN with Inventory System (optional)

Model No.	ML0423WN	ML1626WN	ML3126WN
Temperature range	-18 °C to -25 °C	-18 °C to -25 °C	-18 °C to -25 °C
Forced air circulation	no	no	no
Defrost (automatic / manual)	manual	manual	manual
Reference probe	no	no	no
Electronic controller	yes	yes	yes
Digital display	yes	yes	yes
Temperature alarm visual / audible	yes	yes	yes
Excess and insufficient temp. alarm	yes	yes	yes
Open door alarm	yes	yes	yes
No-volt remote alarm	yes	yes	yes
External finish	painted steel plate	painted steel plate	painted steel plate
Colour	white	white	white
Lock	yes	yes	yes
Door type	solid door	solid door	solid door
Door hinges	changeable	changeable	changeable
Door seal	not changeable	not changeable	not changeable
Internal finish	plastic, white	plastic, white	plastic, white
Interior light	no	no	no
Port hole Ø 16 mm	yes	yes	yes
Spark free	optional	optional	optional
Shelves (number)	1	3	0
Cooling spirals (number)	0	0	5
Shelf positions (number)	1	3	5 (fixed)
Shelf size (mm) (WxDxH)	365x190x6	495x325x6	-
Drawers (number)	0	0	4 and 2 flaps
Wire baskets (number)	0	0	0
Optional drawer system	yes	Yes	yes
Adjustable feet (number & position)	4 feet	4 feet	2 front feet 2 rear rollers
Gross capacity (l)	39	98	232
Net capacity (l)	35	84	212
External width (mm)	480	595	595
External depth (mm)	505	642	634
Door open depth (mm)	905	1130	1185
External height (mm)	565	850	1570
Internal dimensions (mm) (WxDxH)	360x340x395	420x430x640 step (DxH): 165x185	432x430x1300 step (DxH): 190x255
Product weight (kg)	18	37.5	60
Voltage supply	1/N/PE-230 V 50 Hz	1/N/PE-230 V 50 Hz	1/N/PE-230 V 50 Hz
Packed width (mm)	700	700	700
Packed depth (mm)	750	750	750
Packed height (mm)	740	1075	1735
Packed weight (kg)	25	44.5	67
Refrigerant type CFC / HCFC free	R600a	R600a	R600a



ML3126WN



Recommendations for additional storage sets can be found in our separate catalogue. Please ask for a personal fitting option.

