

Racks for ULT / XLT U100



ORGANIZE AND OPTIMIZE THE FREEZER

Quality, efficiency and capacity at the best possible economy, is a vital objective in Nordic Lab's business philosophy. Our storage solutions include racks and boxes, all insuring an optimal handling and market low cost per sample to be stored.

Please find below specifications on the storage solution for the ULT / XLT U100 freezer model.

Rack for 50 mm box



Rack for 100 mm box



Storage capacity

Box height	No. of boxes	No. of racks	Total no. of boxes	Racks without B&D* Product no.
50 mm	74	1	74	N321700
75 mm	45	1	45	N321800
100 mm	33	1	33	N331400

Storage capacity - Rack with boxes and dividers

Box height	Divider	Numbers of samples	Product no.
50 mm	9x9	5994	N322700
50 mm	10x10	7400	N322710
75 mm	9x9	3645	N322800
75 mm	10x10	4500	N322810
100 mm	9x9	2673	N332400
100 mm	10x10	3300	N332410

*Boxes and Dividers