

CYTO-CHAMBERS REF. PV 248 / PV 250

EZ Cytofunnel

For single use with all cyto-centrifuge applications:

- ✓ Optimizes cell performance.
 - ✓ Produces cost-effective, thin-layer preparations
 - ✓ Reduces the risk of exposure to pathogenic samples.
 - ✓ Disposable, to eliminate time-consuming decontamination and clean-up.
 - ✓ Reduces the possibility of cross contamination.
 - ✓ Eliminates the need of a metal slide clip.
- **EZ Single Cytofunnel with white filter card:** Disposable sample chamber with white filter card and lid. Made for faster absorption:
 - Use for sample volumes of 0.1-0.5 ml.
 - Cell deposition area of 6 mm de diameter (28 mm²).
 - White filter card pre-adjusted for high absorption.

