

## Laboratory centrifuge CD-0412-50

- Fix angle rotor (54°) for 6 x 50 ml or 6 x 15 ml
- Automatic self-diagnosis to control all functions
- backlit LCD-display for easy-read out of parameter
- such as Speed, RCF and time
- easy, intuitive operation
- Timer function 0,5 – 99 minutes
- Timer starts after the selected speed is reached
- Quick spin function
- Lid with electric lock and automatic release function to save processing time and prevent sample overheating
- Maintenance free brushless motor
- 2 years warranty



Speed	300-4000 rpm
Increment	100 rpm
max. RCF	1800 g
Timer	0,5-99 minutes, Impulse
continuous operation	yes
Dimension	30x35x20 cm (DxWxH)
Weight	6 kg
max. Capacity	6 x 15 ml or 6 x 50 ml
Rotor	54 ° fix angle
Noise	< 54 dB (A)
Delivery	incl. Rotor
Voltage	100 – 240 V
Frequency	50/60 Hz

