

## Incubation shaker IS-OS20

The new shaking incubator IS-OS 20 combines two typical laboratory operation: in one instrument: shaking and incubation of samples. It is ideal for cell culture, solubility studies, extraction procedures and many other laboratory applications.

The built-in fan guarantees a uniform distribution of the temperature.

IS-OS20 will be supplied with a standard platform equipped with flexible springs that can accommodate many types of flasks, beaker, tubes with different sizes.

We are granting 2 years warranty

### Technical Data:

	IS-OS20
Temperature Range /Resolution	RT+5 -60/0,1°C
Temperature homogeneity at 37°C	±0,5°C
Temperature variation at 37°C	±0,1°C
Timer (min)	99:59/continuous operation
Shaking Speed (rpm)	40....300
Amplitude of platform Movement	20 mm
Capacity	7 x 500ml Erlenmeyer flasks 4 x 1000ml Erlenmexer Flasks
Dimension Platform (WxHxD)	320 x 320 mm
Voltage/Power	230V/500W
External Dimension (WxHxD) mm	500 x 610 x 470
Weight (kg)	40
Adjustable safety temperature	Yes
Material of hood/sickness	Plexiglass/ 8 mm



Phoenix Instrument GmbH – Heinkelstraße 4 – 30827 Garbsen  
[www.phoenix-instrument.de](http://www.phoenix-instrument.de) – [info@phoenix-instrument.de](mailto:info@phoenix-instrument.de)