

## Stirrer RSM-14 HP

### RSM-14HP

- Heating range up to 380 °C
- Work Plate 140 x 140 mm
- Maintenance-free, brushless motor
- Max stirring volume 5 Liter
- Easy-to-read LCD-Display
- digital control of Speed and Temperature
- Connector for external Temperature sensor Pt-1000
- Built-in Temperature control function
- 2 Years warranty



RSM-14 HP with  
PT-1000 and holder



RSM-14 HP

#### Technical Data

Work plate Dimension (mm)	<b>RSM-14 HP</b> 140 x 140 mm
Work Plate Material	Aluminum with porcelain coating
Motor	DC, brushless
Power output [W]	510
Voltage/ Frequency	100-240V/50/60hz
Max. stirring Volume [L H <sub>2</sub> O]	5
Max. length stirring bar [mm]	50
Speed [rpm]	200-1500
Display	LCD
Heating Power [W]	500
Temperature range [°C]	RT+5 - 380
Safety temperature [°C]	420
External Temperature sensor	PT1000
Protection class accor. EN60529	IP 21
Dimension [B x T x H] mm	180 x 320 x 108
Weight	2,2 kg