

Bottle-top dispenser

Choose the right dispenser for your application

		LABmax premium	LABmax eco	LABmax airless	LABmax airless HF	MINISENSOR	TITREX
Suitability	Aqueous solutions	X	X	X	X	X	X
	Crystalline solutions	X	X			X	X
	Solvents	X*	X*	X	X	X*	X*
	HF				X		
Volume setting	Quick-lock knob	X	X			X	
	Adjustment screw			X	X		
Special features	Air-purging	X	X	X	X		X
	180° drip-free system	X	X	X	X		X
	Zero loss of reagent	X	X				X
	360° safety turn	X	X	X	X		X
	Built-in recalibration	X					X
Autoclavable	at 121°C	X	X	X	X	X	X**
Dispensing system	Direct displacement	X	X			X	X
	Contactless piston			X	X		
Certifications	Individual certificate	X	X	X	X	X	X
	DE-M marked	X	X	X	X	X	X

*For detailed information regarding the chemical resistance of the dispenser, please see page 16.

**Only valve block autoclavable.

