

The economic electric pipetting solution

WITOPED eco is an economically priced electrical pipette aid, suitable for all pipettes up to 100 ml.

Nosepiece

The nosepiece is suitable for all kinds of pipettes. In case of contamination, the nosepiece can be easily removed in order to exchange the filter and the silicone adapter.

Dispenser buttons

The WITOPED eco has two buttons, one for aspirating the liquid (up), one for releasing the liquid (down).


Digital display and speed setting

The LCD display indicates battery status and set speed. The pipetting speed can be adjusted by pushing the "+/-" buttons.



Motor

The powerful motor ensures fast and reliable pipetting.





Pipetting aid

Your advantage at a glance:

- suitable for all pipettes up to 100 ml
- most affordable price in its class
- digital display indicating battery status and selected speed
- efficient battery for a long runtime, low battery alert
- fine adjustment of pipetting speed
- tight fit of any pipette into the pipette socket
- handy and light weighted for fatigue-free pipetting
- simple but powerful motor for fast and reliable pipetting
- gravity drain mode for draining out without using any motor
- nosepiece with replaceable filter and silicone adapter
- silicone adapter autoclavable at 121°C
- including Li-ion 3,7V 1100mAh battery
- CE label

WITOPED eco

Type	Order number
220 V	5 381 000

Scope of supply


WITOPED eco, battery charger, replacement filter, table stand, instruction manual

Replacement one-way filter



Packaging unit	Order number
5	5 380 006 E

Replacement silicone adapter



Order number
5 382 010