


“Witeg” **WiseTis® HG-15D** Digital Homogenizer with Remote Controller, 1~2,500ml
With Microprocessor Digital-Control and LCD display of Speed & Power, 2,000~2,7000rpm, "CE-Certificate"



HG-15D-Set

Applications & Features

- Designed for Liquid and Tissue Samples
- Useful for All Dispersing and Homogenizing Applications.
- High-speed / High Quality / Low Noise :
Max. 27,000 RPM
- Remote Digital controller :
Fully Digital control using Microprocessor
- Speed & Power Control optimized for the Sample
- LCD Display of Speed(RPM) and Power Bar Graph
- High Efficiency Dispersing Tool made of Stainless Steel and PTFE
- Over-Load & Over Heat Protection Functions
- Low Cost & High Efficiency

Specification

Model		HG - 15D
Dimension (w×d×h mm)		60 × 70 × 230
Working Range (ml)		1 ~ 2,500
Speed (rpm)		2,000 ~ 27,000
Speed Control Resolution		10 rpm
Motor Rating (Input / Output)		300 / 160 W
Controller		Digital Feedback Controller (Remote type)
Display		Digital LCD Display (Timer, RPM and Power bar - graph)
Usable Tool		HT1008, HT1010, HT1018 & HT1025
Material	Body	Powder Coated Steel
	Tool	Stainless Steel (#304)
Safety Device		Over Load & Over Heat Protector
Weight (kg)		2.5
Power Supply *		AC 220V, 50/60Hz

* Others Specifications are available upon customer's request.

Ordering Information

Homogenizer HG-15D-Unit, with Remote Controller, Digital

DH.WHG02018 Homogenizer HG-15D-Unit, 1~2500ml
with Remote Controller, without Dispersing Tool

Homogenizer Package, HG-15D-Set, Complete with Remote Digital Controller
Dispersing Tool, Clamp and Stand.

DH.WHG02019 Homogenizer Package, HG-15D-Set
with Remote Controller, HT1018 Dispersing Tool,
CL200 clamp and ST110 Standard Plate Stand

DH.WHG02020 Homogenizer Package, HG-15D-Set
with Remote Controller, HT1025, Dispersing Tool,
CL200 clamp, ST110 Standard Plate Stand

Dispersing Tools

DH.WHG02108 Homogenizer Tool, HT1008
Rotor ϕ 6.1, Stator ϕ 8, 1~50ml

DH.WHG02110 Homogenizer Tool, HT1010
Rotor ϕ 7.5, Stator ϕ 10, 1~100ml

DH.WHG02118 Homogenizer Tool, HT1018
Rotor ϕ 12.7, Stator ϕ 18, 10~1000ml

DH.WHG02125 Homogenizer Tool, HT1025
Rotor ϕ 17, Stator ϕ 25, 50~2500ml

Additional Function Accessories

※ DAIHAN WiseStir®, WiseTis® Accessories

DH.WOS502020 CL200, Clamp for Stirrer or Homogenizer
Cast Aluminum, Grip Capacity ϕ 23mm
for Stand/Rod

DH.WOS501510 ST110, Standard Plate Stand for Stirrer and
Homogenizer, ϕ 16 × h800 mm Rod
20 × 35 cm Plate, with Silicon Pad

DH.WOS501515 ST120, Wide Plate Stand for Stirrer and
Homogenizer, ϕ 23 × h1000 mm Rod
30 × 40 cm Plate, with Silicon Pad

※ See page 31 for the specifications of Dispersing tools and Accessories.