

# Laminar-flow cabinet BIO130 CYTO

🏠 [Laminar airflow Safety class II > Laminar-flow cabinet BIO130 CYTO](#)



## Description

- vertical laminar air flow
- 3 HEPA filters with efficiency of 99.995% for MPPS
- microprocessor control system with touchscreen display in Polish/English or any other language
- automatic filter compensation ensures constant and safe speed of air flow in the cabinet
- hour timer of the cabinet and UV lamp operation
- adjusted air flow
- clear and easy-to-use control panel with optical and acoustic signals
- silent fans, noise level below 56 dB
- front window made of tempered double glass without frame, inclined at an angle of 8°, electrically raised and lowered, with a possibility to shut it completely – closing the workspace
- split worktop (4 autoclavable boards), made of stainless steel
- internal workspace made of stainless steel
- permanently fixed UV lamp in the upper part of the workspace, protected from accidental switching on during operation
- white LED lightning system min. 1200 lux
- workspace sides with window openings
- two power sockets inside workspace
- removable armrest
- frame with lockable castors
- power supply 230V/60 Hz

## Gallery



## Warranty and certificates

Product and manufacturer certificates:

- TUV certificate for compliance with the standard EN 12469
- CE declaration of the manufacturer
- ISO 9001



Protocol from tests in laboratory is issued and enclosed for each cabinet.

Warranty for the laminar airflow unit: 36-month warranty.





<b>Nominal width</b> in mm	1340	
<b>External dimensions</b> width/depth/height in mm	1333 / 790 / 2069	
<b>Internal dimensions</b> width/depth/height in mm	1266 / 644 / 710	
<b>Area of working space</b> in m <sup>2</sup>	0,81	
<b>Number of perforated plates in worktop</b> in pcs.	4	
<b>Maximum height of window panel lifting during work</b> in mm	200	
<b>Maximum height of front window panel lifting</b> in mm	440	
<b>Air flow</b> in m/s	inlet	0,53
	working space	0,35
<b>Air flow</b> in m <sup>3</sup> /h	inlet	376
	operation 63%	761
	outlet 37%	376
<b>Number of filters HEPA or ULPA</b> in pcs.	preliminary	1
	main	1
	outlet	1
<b>Filter efficiency</b>	> 99,995 % for MPPS - H14	
<b>Loudness level</b> in dB	< 56	
<b>Illuminance</b> in lux	> 1200	
<b>Maximum power consumption</b> in W	400	
<b>Power of UV bactericidal lamp</b> in W	18	
<b>Number of electric sockets inside the chamber</b> standard, in pcs.	2	
<b>Control panel</b> with touchscreen	menu in Polish, English or other language	
<b>Gas valve</b> (option)	mechanical, with 3/8" connector	
<b>Chamber construction</b>	Floor, working plates, sides and back of the working space	Steel sheet AISI 304(1.403). Acid-proof, galvanized steel sheet
	Body	Steel 1203, thickness 1,5 mm, powder coated RAL 7035
	Front window panel	Safety glass, thickness 6 mm, inclined at an angle of 8°
<b>Equipment weight</b> in kg	280	

In the interest of continued product development, ALPINA reserves the right to make improvements in this leaflet and the products it describes at any time, without notices or obligation. Colors of products may not be perfectly accurate due to variations caused by photograph, printing method or monitor settings.

