

Placa fría para histología | **Cooling plate for histology**

Marca|**Trademark: Nahita**

La Placa Fría Nahita permite el enfriamiento y solidificación de los bloques de parafina líquida obtenidos previamente en la unidad dispensadora de parafina, completando de esta manera el proceso de preparación de muestras de tejido para su posterior corte en el microtomo.

Principales características:

Placa de enfriamiento de gran eficacia, con borde de contención metálico.

- Espacio en placa para un máximo de 57 bloques
- Enfriamiento en continuo
- Libre de CFC

Liquid paraffin blocks obtained from the dispensing unit can be cool thanks to Nahita cold plate, completing the preparation process of tissue samples for their subsequent cut with the microtome.

Main features:

- Cold plate with high efficacy, metallic edge.
- Maximum space for 57 blocks.
- Constant cooling
- Free of CFC



Referencia / <b>Reference</b>	50550002
Capacidad / <b>Capacity</b>	Hasta 57 bloques Max. 57 blocks
Temperatura / <b>Temperature</b>	Hasta -20°C±
Dimensiones (placa) / <b>Dimensions (plate)</b>	310x298 mm
Dimensiones (equipo) / <b>Dimensions (equipment)</b>	570x345x385 mm
Alimentación / <b>Supply</b>	AC 220±10% (50Hz)
Peso / <b>Weight</b>	23 Kg