



thermoblock  
nahita  
MINIT

thermoblock, Mini T

## **thermoblock** **Mini T**

New Mini “T” thermoblock easy to use and designed with special compact dimensions that fit in the palm of your hand. Available several interchangeable aluminum blocks and precise control of temperature that gives versatility for use in biochemistry research, medical, chemical or veterinary laboratories being very useful in molecular biology and genetics (DNA samples, extraction of lipids or other cellular components, denaturalization of proteins, etc.).



thermoblock, Mini T

## Thermoblock features

- [01] LCD display easy to set up and use
- [02] Simultaneous display of sent and actual time and temperature.
- [03] Temperature deflection adjusting meet with different demands.
- [04] Beep-signal at the end of the program.
- [05] Mini-dimensions: very compact design that fits in the palm of your hand.
- [06] Several aluminum blocks (not included) are available covering a wide range of tubes sizes (from 0.2 ml to 50 ml) and other type of recipients.



## Thermoblock back of the equipment



- [01] Switch on/off
- [02] Current connexion

thermoblock, Mini T

## Thermoblock technical specifications



Code	50603000
Model	Mini T
Temperature range	RT+5 °C – 80 °C
Time	1 – 999 s / 1 – 999 Min
Temperature control accuracy	± 0.5 °C
Temperature accuracy	0,1 °C
Heating time (de 20°C a max. Temp.)	≤12 min
Temp. control accuracy at 40°C	0,3 °C
Dimensions (L x W x H)	110 x 150 x 80 mm
Weight	0,5 Kg

## Thermoblock acesories (not included)

Co	Number of tubes	Capacity
90603001	40	0,2 ml
90603002	24	0,5 ml
90603003	15	1,5 ml
90603004	15	2 ml
90603005	8 cells	12,5 x 12,5 mm
90603006	4	15 ml
90603007	2	50 ml

[01]



[02]



[04]



[03]

[05]



[06]



[01] 90603001

[02] 90603002

[03] 90603003

[04] 90603005

[05] 90603006

[06] 90603007

AUXILAB, S.L.

Polígono Morea Norte Calle D Nº 6

31191 Beriain - Navarra - España

T: +34 948 310 513

comercial@auxilab.es

www.auxilab.es

ROGO-SAMPAIC

23 Rue Ampère, ZI de Vilemilan

91325 Wissous - Paris - France

T: +33 1 69 53 67 67

rogosampaic@rogosampaic.com

www.rogosampaic.com



AUXILAB S.L.

Material de laboratorio  
Laboratory supplies