

Baño termostático par cortes de parafina | **Thermostatic bath for parffin sections**

Marca|Trademark:Nahita

Novedoso equipo en diseño y materiales, de alta estabilidad funcional y fácil manejo, que conforma parte del instrumental empleado en Histología y Patología para el tratamiento de tejidos. Características principales:

- Amplio panel de control de fácil comprensión y manejo
- Cubeta termostaticada para la solución de montaje que, a la temperatura adecuada, permite la separación y posterior aplicación en el portaobjetos de la tira de cortes obtenida mediante el microtomo.
- Superficie en negro de la cubeta que permite una buena observación del corte.
- Amplio borde térmico para un excelente secado de la preparación.

Equipment with new design and materials, high functional stability and easy to use, it is used in Histology and Pathology for tissues treatment.

Main features:

- Wide panel control, which is easy to understand and use.
- At the appropriate temperature the thermostatic chamber permits the separation and subsequent application on the microscope slide of the strip obtained with the microtome.
- Black surface of the chamber for a good observation of the section.
- Wide thermal limit to get an excellent drying of the preparation.



Referencia / Code	50605001
Capacidad / Capacity	~2L
Rango de Temperatura / Temperature range	Ambiente hasta 90°C / Room temperature to 90°C
Precisión / Precision	-3 / +5 °C
Dimensiones (cubeta) / Dimensions (chamber)	240x180x54 mm
Dimensiones (equipo) / Dimensions (equipment)	340x350x105 mm
Alimentación / Supply	AC 220±10% (50Hz)
Consumo / Consume	550 w
Fusibles / Fuses	6A
Peso / Weight	4.5 kg