

BOECO TUBE ROLLER TR-6

the microprocessor controlled BOECO TR-6 provides gentle rocking and rolling motion. The user can define programs (speed & time) and its Pulse mode provides a momentary pause for through mixing (Pulse can be set for a range of 30 sec to 90 sec. time) The rollers are fully autoclavable and easy to remove which could allow to use bigger tubes.

Specification:


Housing:	Powder coated sheet metal
Motor:	Brushless DC motor
Speed range:	10 - 80 rpm
Speed display:	LED
Timer:	1 minute to 99 min, / non-stop
Roller:	6 Rolls / 5 positions (Nylon - autoclavable)
Roller Length:	327 mm
Overall size:	496 x 292 x 107 mm
Weight:	approx. 6 kg

Code

BOE 8025000

Description

Tube roller TR-6, AC adapter 100-240V, 50/60Hz



TR-6