

# Table top Centrifuge LACE32



## λ Application

Our clinical centrifuge is ideal for the separation of serum, plasma, urea, blood samples and other routine applications in hospital and research laboratories.










## λ Features

1. Brush-less frequency motor drive, maintenance is free
2. Table top low speed multipurpose centrifuge
3. Melt housing, 304# stainless steel centrifugal chamber, safe and reliable.
4. Microprocessor control and digital display run parameters, like speed, RCF, time.

## TECHNICAL SPECIFICATION

<b>Max.RPM/RCF</b>	<b>Angle rotor</b>	<b>5,000/3860×g</b>	
	<b>Swing out rotor</b>	<b>5,000/4730×g</b>	
<b>Max.capacity</b>	<b>Angle rotor</b>	<b>12×100ml</b>	
	<b>Swing out rotor</b>	<b>4×250ml</b>	
<b>Timer</b>	<b>1~99min</b>	<b>RPM/RCF conversion</b>	<b>Yes</b>
<b>Noise level (dB)</b>	<b>≤55</b>	<b>Acc/Dec (sec)</b>	<b>10 steps</b>
<b>Digital display</b>	<b>Yes</b>	<b>Safety lid lock system</b>	<b>Safe lid lock</b>
<b>Power supply</b>	<b>220V 60Hz</b>	<b>Temperature range</b>	<b>Normal system</b>
<b>Product Dimension (W×D×Hmm)</b>	<b>558×400×300mm</b>	<b>Package Dimension (W×D×HMM)</b>	<b>635×460×360mm</b>
<b>Weight (kg)</b>	<b>N.W:35kg</b>	<b>G.W:38kg (carton box)</b>	
<b>Certificates</b>	<b>CE,ISO &amp; Calibration report are available</b>		

## OPTIONAL ADAPTERS AND ROTORS FOR LACE32

	<b>Rotor No.</b>	<b>Volume (ml)</b>	<b>Max.Speed (RPM)</b>	<b>Max.RCF (×g)</b>
	30569	4×250ml	4000	2880
	31277	32×10ml	4000	2880
	31277	24×15ml	4000	2950
	31377	48 ×7/5ml	4000	2760/2480
	30641	12×50ml	5000	3860
	30614	6×100ml	5000	3130
	30642	12 ×100ml	4000	2970
	30639	24×15ml	5000	3500
	31376	2× 3×48well	4000	2300