

ACTIVE COAL FILTERS SELECTION GUIDE

Flammable cabinets

Corrosive cabinets

Toxic cabinets and multirisk








Pesticide range

Fume hoods, filtering cabinets - Ventilation

Retention and Cans










Anti-fire equipments and Files cabinets

Showers and first-aid equipments

	Kind of filter	Filter for vapors	Reference from 2011	For	Reference before 2011
Filtration boxes					
 H 150 x Ø 150 (mm)		Organic	CDFORG	Range V - Boxes: CDF-A, CDF-P, CVAT, CMAT, CTAT	CDF-F20
		Corrosive	CDFCORG		CDF-F35
		Polyvalent ⁽¹⁾			CDF-F15
		Formaldehyde	CDFFOR		-
 H 70 x L 390 x P 390 (mm)		Organic	ORG50	Range H - Boxes: H50	ORG50
		Polyvalent ⁽¹⁾	CORG51		CORG51
		Formaldehyde	FOR50		-
Filtering ventilation cabinets					
 H 150 x Ø 150 (mm)		Organic	CDFORG	Range 12 Cabinets: AF10	CDF-F20
		Corrosive	CDFCORG		CDF-F35
		Polyvalent ⁽¹⁾			CDF-F15
		Formaldehyde	CDFFOR		-
 H 70 x L 390 x P 390 (mm)		Organic	ORG50	Range 12 - Cabinets: AF1, AF2 (before February 2002)	CDF720
		Corrosive	CORG51		CFR735
		Polyvalent ⁽¹⁾			CFR - CFR70 - CFR715
		Formaldehyde	FOR50		-
		 H 20 x L 390 x P 390 (mm)			Polyvalent ⁽¹⁾
 H 200 x L 100 x P 400 (mm)		Organic	ORGFC	Range 12 - Cabinets: AF1, AF2, AF2-2, AF4-2, AF4-3, AF4-4 (before July 2008)	FC20
		Corrosive	CORGFC		FC35
		Polyvalent ⁽¹⁾			FC15
		Formaldehyde	FORFC		-
 H 70 x L 730 x P 380 (mm)		Organic	ORG300	Range 12.X - Cabinets: AF1X, AF2,X AF2-2X, AF4-2X, AF4-3X, AF4-4X (after July 2008)	ORG300
		Corrosive	CORG301		COR300
		Polyvalent ⁽¹⁾			CORG301
		Formaldehyde	FOR300		-

⁽¹⁾ organic + corrosive. CORG polyvalent active coal filters can adsorb the most common chemical products used: organic and inorganic acids, organic compounds or solvents.

ACTIVE COAL FILTERS SELECTION GUIDE

	Kind of filter	Filter for vapors	Reference from 2011	For	Reference before 2011		
Fume Hoods							
	 H 70 x W 390 x D 390 (mm)	Organic	ORG50	Range H Fume hoods : H60	-		
		Polyvalent ⁽¹⁾	CORG51		-		
		Formaldehyde	FOR50		-		
		Ammonia	AMM50		-		
	 H 70 x W 585 x D 380 (mm)	Organic	ORG200	Range H Fume hoods : H200, H201 H210, H211 PMSF200	ORG200		
		Corrosive	CORG201		COR200		
		Polyvalent ⁽¹⁾	CORG201		CORG201		
		Formaldehyde	FOR200		-		
		Ammonia	AMM200		-		
		Organic	ORG300 ORG400 ORG1300 ORG2300		Range H Fume hoods : H300, H301, H310, H311, H400, H401, H410, H411, H1300, H1301, H1310, H1311, H2300, H2301, H2310, H2311, PMSF400, PMSF1300, PMSF2300	ORG300 ORG400 ORG1300 ORG2300	
		Corrosive	CORG301 CORG401 CORG1301 CORG2301			COR300 COR400 COR1300 COR2300	
Polyvalent ⁽¹⁾	CORG301 CORG401 CORG1301 CORG2301	CORG301 CORG401 CORG1301 CORG2301					
Formaldehyde	FOR300 FOR400 FOR1300 FOR2300	- - - -					
Ammonia	AMM300 AMM400 AMM1300 AMM2300	- - - -					
Organic	ORG500 ORG1500 ORG2500 ORG3500	Range H Fume hoods : H500, H501, H510, H511, H1500, H1501, H1510, H1511, H2500, H2501, H2510, H2511, H3500, H3501, H3510, H3511, PMSF1500	ORG500 ORG1500 ORG2500 ORG3500				
Corrosive	CORG501 CORG1501 CORG2501 CORG3501		COR500 COR1500 COR2500 COR3500				
Polyvalent ⁽¹⁾	CORG501 CORG1501 CORG2501 CORG3501		CORG501 CORG1501 CORG2501 CORG3501				
Formaldehyde	FOR500 FOR1500 FOR2500 FOR3500		- - - -				
Ammonia	AMM500 AMM1500 AMM2500 AMM3500		- - - -				
	  (2x) H 70 x W 585 x D 380 (mm)		Organic	ORG500 ORG1500 ORG2500 ORG3500	Range H Fume hoods : H500, H501, H510, H511, H1500, H1501, H1510, H1511, H2500, H2501, H2510, H2511, H3500, H3501, H3510, H3511, PMSF1500	ORG500 ORG1500 ORG2500 ORG3500	
			Corrosive	CORG501 CORG1501 CORG2501 CORG3501		COR500 COR1500 COR2500 COR3500	
		Polyvalent ⁽¹⁾	CORG501 CORG1501 CORG2501 CORG3501	CORG501 CORG1501 CORG2501 CORG3501			
		Formaldehyde	FOR500 FOR1500 FOR2500 FOR3500	- - - -			
		Ammonia	AMM500 AMM1500 AMM2500 AMM3500	- - - -			
		 H 100 x W 565 x D 405 (mm)	Organic	ORG203		Compliant with MF3-AS	ORG203
			Corrosive	CORG203		Compliant with MF3-BE	CORG203
			Polyvalent ⁽¹⁾	CORG203		Compliant with MF3-F	-
			Formaldehyde	FOR203		Compliant with MF3-K	-
			Ammonia	AMM203		Compliant with MF4-AS	ORG303
Organic	ORG303		Compliant with MF4-BE	CORG303			
 H 100 x W 755 x D 405 (mm)	Corrosive	CORG303	Compliant with MF4-F	-			
	Polyvalent ⁽¹⁾	CORG303	Compliant with MF4-K	-			
	Formaldehyde	FOR303					
	Ammonia	AMM303					
		Specific, contact-us.	Don't forget, we are able to provide specialised filters for specific applications for fume hoods, filtering ventilation cabinets, filtration boxes...in standard manufacturing and in specific manufacturing, contact-us for more information.				

⁽¹⁾ organic + corrosive. CORG polyvalent active coal filters can adsorb the most common chemical products used: organic and inorganic acids, organic compounds or solvents.

Flammable cabinets

Corrosive cabinets

Toxic cabinets and multirisk

Pesticide range

Fume hoods, filtering cabinets - Ventilation

Retention and Cans

Anti-fire equipments and Files cabinets

Showers and first-aid equipments

The manufacturer keeps the right to modify technical characteristics without notice. 12/2013