

COOL'S EB04.003

Isothermal bag for the transportation of samples



High capacity isothermal bag that allows the transport, in an organized way and under optimal conditions, of 300 diagnostic samples and 36 urine samples, from the extraction center to the laboratory.

The bag is certified according to UN 3373, that enables the transport of biological substances class B.

Optionally EB04.003 Cool's bag can be delivered with 6 or 12 pieces of reusable cooling gel packs EB09.007 (25,5 x 10,5 cm)

*** Final production may vary from the one shown in the pictures.**

Features

- This bag consists of:
 - Main isothermal compartment
 - 4 detachable isothermal compartments:
 - 1 special compartment for urine containers, with removable dividers. It can hold up to 24 containers. Size: 42x26.5x9.5cm.
 - 3 compartments for test tube racks. Each one has a capacity for 2 racks (100 test tubes), with a size of 26.5

[wEBsites](#)

All detachable compartments have a handle, identifier and

Technical data

- **Measures:** 44 x 29 x 39 cm
- **Capacity:** 49,76 L
- **Weight:** 5.30 kg
- **Material:** Tarpaulin. Water resistant.
- Washable
- **Colour:** Blue
- **Maximum recommended contents weight:** 10 kg

Key features

- Upper transparent pocket for clinical history and documents.
- Transparent zippered front pocket
- Handle, shoulder strap and trolley band.
- Each detachable compartment has a special flap to place the cold gels.



tee-uu

social media

