

# Universal oven "Paraffin Oven"

## Model UN55pa

### Standard equipment

#### Ventilation and Control

- natural convection
- adaptive multifunctional digital PIDmicroprocessor controller with 2 highdefinition TFT-colour displays
- self-diagnostics for fault analysis
- 2 Pt100 sensors Class A in 4-wire-circuit, mutually monitoring and taking over the performance at the same temperature value
- digital timer adjustable from 1 minute to 99 days
- parameters adjustable on ControlCOCKPIT: temperature (Celsius or Fahrenheit), programme time, time zones, summertime/wintertime
- SetpointWAIT function the process time does not start until the set temperature is reached
- HeatBALANCE function for adapting the distribution of the heating performance of the upper and lower heating circuit from -50% to +50%
- Setting of languages on ControlCOCKPIT: German/English/Spanish/French
- resolution of display for setpoint and actual temperature values 0.1 ℃
- calibration facility directly at ControlCOCKPIT for three freely selectable temperature values
- programme stored in case of power failure
- incl. works calibration certificate for +80 ℃

#### Multiple Overtemperature Protection

- with audible and visual alarm
- overtemperature monitor TWW, protection class 3.1 or adjustable temperature limiter TWB, protection class 2, selectable on display
- additionally integrated over- and undertemperature monitor "ASF", automatically following the setpoint value at a preset tolerance range; alarm in case of overor undertemperature, heating is switched off in case of overtemperature
- mechanical temperature limiter TB, protection class 1 according to DIN 12 880 to switch off the heating approx. 20 °C above nominal temperature

#### **Textured Stainless Steel Casing**

- w x h x d: 585 x 787 x 514 mm
- fully insulated stainless steel door with 2-point locking (compression door lock)
- rear zinc-plated steel

#### **Interior - Heating Concept**

- w x h x d: 400 x 400 x 330 mm, 53 l
- easy-to-clean interior,made of stainless steel,reinforced by deep drawn ribbing with integrated and protected large-area heating on four sides
- nearly gas-tight interior
- 1 stainless steel grid
- Max. loading of chamber: 80 kg

#### **Temperature Range**

• min. 5 ℃ above ambient up to +80 ℃

#### Voltage / Power Rating

- 230 V (+/- 10%), 50/60 Hz / approx. 2.000 W
- 115 V (+/- 10%), 50/60 Hz / approx. 1.700 W

#### **Packing Data**

- net weight approx. 57 kg
- gross weight carton approx. 69 kg
- dimensions approx.: carton w x h x d: 70 x 104 x 73 cm
- the appliances must be transported upright

#### **Customs Tariff Number**

• 8419 8998

#### **Country of Origin**

• Federal Republic of Germany

#### WEEE-Reg.-No.

• DE 66812464