

<b>Biocen 22 / R</b>		<b>RT 222</b> 		<b>RT 223</b> 		<b>RT 224</b> 	
<b>ROTOR</b>		<b>Angle fixed 45°</b>		<b>Horizontal</b>		<b>Angle fixed 45°</b>	
Capacity		24 de 1,5x2 ml.		8x15 ml. conical		32 de 0,2 ml.	
R.P.M.		18100		8000		18100	
Radius mm.		82		91		55	
R-C.F.		30034		6511		20145	
Min.temp.max speed.(°C)		0		-2		0	
		<b>ADAPTERS</b>		<b>ADAPTERS</b>		<b>ADAPTERS</b>	
<b>Tubes capacity</b>	<b>Dim (mm)</b>	<b>TUBES</b>	<b>REF.</b>	<b>TUBES</b>	<b>REF.</b>	<b>TUBES</b>	<b>REF.</b>
15 ml.	16x100	-	-	8	-	-	-
10 ml.	13x100	-	-	8	RE 371	-	-
5 ml. blood test	13x82	-	-	8	RE 377	-	-
1,5x75 mm. capillarie	1,5x75	-	-	-	-	-	-
microtubes 1,5-2 ml.	11x39	24	-	-	-	-	-
microtubes 0,5-0,6 ml.	8x30	24	RE 305	-	-	-	-
microtubes 0,2-0,4 ml.	6x45	24	RE 304	-	-	32 de 0,2	-