

## High-Speed Mini Centrifuge

UEC 15



# UNIEQUIP



- Max Speed of 15.000 rpm > RCF 15345
- Rotor 18 x 1,5/ 2 ml with aerosol tight lid
- Reduction adapters for 0,2/0,4/0,5 ml
- PCR Rotor 2 x 8 > 0,2 µl optional
- Imbalance detection autom. out/off
- Timer function 0-999 minutes& hours
- Big black & white LCD display
- USB port for documentation
- Small footprint , low noise

Dimensions and specifications are subject to change without notice for further improvements on design.  
Dimensions may vary dependent on different accessories be used.

## UNIEQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11  
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04  
e-mail: [info@uniequip.de](mailto:info@uniequip.de) · Internet: [www.uniequip.com](http://www.uniequip.com)