

“witeg” **Wisemix® BML** Programmable Ball Mill, 2-/6- Pots
with Digital Feedback Control System, Digital LCD Display, 50~600 rpm, ±5 rpm, 10 Steps 6 Programs



2-Places Programmable Ball Mill “BML-2”

- Ideal for Grinding/Mixing : Crystalline Chemicals, Inks, Paints, Pigments, Plants, Racks and Soils
- Microprocessor, Programmable, Variable Speed : 50~600 rpm
- Digital Feedback Control and LCD Display : Speed (RPM), Power Bar-Graph and Timer Status
- CE Certified

- * Antistatic Rubber Roller : Keeps Jars / Pots Centered
- * Digital LCD with Back-Light Function
- * Over Load Protection
- * Ergonomic and Compact Design
- * NBR Roller with Pot Guide Ring
- * Programmable Operation : 10 Steps 6 Programs
- * Storage Function (RPM, Timer)
- * Soundproof Cabinet for Safety and Minimum Noise (Optional)



6-Places Programmable Ball Mill “BML-6”

Specification

Model		2-places“BML-2”	6-places“BML-6”
Capacity		up to 2 pots	up to 6 pots
Motor		DC 40W Geared Motor (5 : 1 ratio)	DC 100W Geared Motor (5 : 1 ratio)
Maximum Load		40kg	120kg
Speed Range		50 ~ 600 rpm	
Speed Accuracy		± 5 rpm	
Timer		99 hr 59 min (Programmable)	
Display		Digital LCD Display with Back-Light Function	
Controller		Programmable Digital Feedback Controller with Jog-Shuttle Switch (Turn+Push)	
Program Function		10 Steps, 6 Programs	
Safety Device		Over Load Protector for Motor	
Dimension (w×d×h)	Roller (mm)	50×300	70×800
	External (mm)	620×300×515	1,200×500×795
Packing Size (w×d×h) (mm)		760×460×700	1,300×620×1,010
& Gross Weight		52 kg	115 kg
Power Consumption		40 W	100 W
Power plug type		Hybrid E+F	
Power Supply*		AC 230V, 1N~, 50/60Hz	

* Other Specifications are available upon Customer's Request.

Ordering Information

For 2 Pots

DH.WBM05002 Ball Mill, 2-places “BML-2”, 230V, with Cord/Plug

For 6 Pots

DH.WBM05006 Ball Mill, 6-places “BML-6”, 230V, with Cord/Plug

Optional Soundproof Cabinets

* Designed for Minimizing Noise, and Safety Work

DH.WBM050102 Soundproof Cabinet, “BMC2”, for 2-places“BML-2”, w1.020×d570×h1.120mm

DH.WBM050106 Soundproof Cabinet, “BMC6”, for 6-places“BML-6” w1.600×d770×h1.400mm



Optional Soundproof Cabinet “BMC-6” for 6-Places Programmable Ball Mill “BML-6”